

A Validated High-Performance Thin-Layer Chromatography Method for Monitoring Glucose and Malto Oligosaccharides with Multiple Degree of Polymerization During Bioethanol Production from Corn Biomass.

McKenzie Britt and Wilmer Perera

CAMAG Scientific, Inc. Wilmer.Perera@camag.com

Introduction

Fossil fuel

Is a hydrocarbon-containing material formed underground from the remains of dead plants and animals.



Coal



Petroleum



Natural gas

Introduction

Advantages

- Fossil fuels generate a large amount of electricity at a single location.
- Can be found easily and are cost-effective.
- Transportation of oil and gas can be done easily through pipelines.

Disadvantages

- Fossil fuels are the main driver of global warming, then contribute to climate change.
- Fossil fuels are non-renewable sources of energy.

Renewable sources of energy

- Solar energy
- Geothermal energy
- Wind energy
- Hydropower from flowing water
- Biomass from plants



Biomass from plants



 Corn has been widely used to produce ethanol due to the large amount of carbohydrates (specifically starch) produced in the endosperm of the kernel.

Amylose







Biomass from plants



 Bioethanol production is a process of sugar fermentation to yield ethanol

• Corn requires from **60 to 100** days to reach harvest depending on the variety and warm weather.

Ethanol production from corn



DRY MILL ETHANOL PROCESS

HPTLC Association

The International Association for the Advancement of High-Performance Thin-Layer Chromatography

Monitoring sugars by HPLC

- The fermentation process of most bioethanol plants in U.S is monitored by HPLC.
- The HPLC methods range from 21 to 42 min and the detection is usually performed using a refractive index detector.
- The HPLC method can simultaneously detect ethanol, organic acids, glucose, maltose and maltotriose.
- Maltotetraose and sugars with the highest degree of polymerization are not resolved by HPLC.



The aim of this work is to develop an HPTLC method to quantify from glucose up to DP₈ and improve the turnaround of the process.

Glucose	6-48 hours= 100x dilution		
	54-60 hours= 5x dilution		
Maltose	No dilution		
Maltotriose	No dilution		
554	6, 36-60 = 2x dilution		
DP4	12-30 hours = 5x dilution		
005	2x dilution		
DP5	6 &12 hours = 5x dilution		
	no dilution		
DP6	6 hours 5x dilution		
DP7	No dilution		
DP8	No dilution		



HPTLC instruments and step process



HPTLC validation

Stability

- Sample in solution (3h)
- Sample on plate (3h)
- 2D-Separation
- Evaluate chromatogram immediately and during 1 hour after derivatization or final drying.



Stability

Stability

Stability of the sample in solution

and on the plate (3h)



2D-Separation



ng

Stability

Evaluation of the chromatogram immediately and during 1 hour after derivatization.





Repeatability

A586-031921-01



A586-031921-02



A586-032021-01



Sugars	A586-031921-01	A586-031921-02	A586-032021-01	$\Delta R_{\rm F}$
Glc	0.64	0.65	0.65	0.01
DP2	0.56	0.57	0.57	0.01
DP3	0.46	0.47	0.47	0.01
DP4	0.36	0.37	0.37	0.01
DP5	0.27	0.28	0.28	0.01
DP6	0.19	0.20	0.20	0.01
DP7	0.16	0.16	0.17	0.01
DP8	0.14	0.13	0.14	0.01

Repeatability









Recovery for glucose





Recovery (%) = <u>S – U</u> × 100

 $\frac{1}{16} \frac{1}{R^2 D} \frac{1}{R_1} \frac{1}{R_2} \frac{1}{R_1} \frac{1}{R_2} \frac{1}{R_1} \frac{1}{R_1} \frac{1}{R_2} \frac{1$

Α

^	Sample 'Glu-0407	21-01-01'	23.64 µg/ml	CV=5.47 %	(3 applications)
	 Volume: 3.0 μl 		23.64 µg/ml	CV=5.47 %	(3 replicas)
	Track 7	RF 0.640	22.91 µg/ml	68.72 ng	
	Track 8	RF 0.642	22.89 µg/ml	68.67 ng	
	Track 9	Rr 0.642	25.14 µg/ml	75.41 ng	
^	Sample 'Glu-0407	21-01-03	50.17 µg/ml	CV=1.23 %	(3 applications)
	 Volume: 3.0 µl 		50.17 µg/ml	CV=1.23 %	(3 replicas)
	Track 10	RF 0.644	49.49 µg/ml	148.5 ng	
	Track 11	Rr 0.644	50.70 µg/ml	152.1 ng	
	Track 12	RF 0.646	50.33 µg/ml	151.0 ng	
^	Sample 'Glu-0407	21-01-04	73.08 µg/ml	CV=0.52 %	(3 applications)
	 Volume: 3.0 μl 		73.08 µg/ml	CV=0.52 %	(3 replicas)
	Track 13				
	Theorem 10	Rr 0.646	72.64 µg/ml	217.9 ng	
	Track 14	Rr 0.646 RF 0.648	72.64 µg/ml 73.35 µg/ml	217.9 ng 220.0 ng	
	Track 14 Track 15	Rr 0.646 Rr 0.648 Rr 0.648	72.64 µg/ml 73.35 µg/ml 73.25 µg/ml	217.9 ng 220.0 ng 219.8 ng	
^	Track 14 Track 15 Sample 'Glu-0413	Rr 0.646 Rr 0.648 Rr 0.648 21-01-05'	72.64 µg/ml 73.35 µg/ml 73.25 µg/ml 95.19 µg/ml	217.9 ng 220.0 ng 219.8 ng CV=2.20 %	(3 applications)
^	Track 14 Track 15 Sample 'Glu-0413 ^ Volume: 3.0 µl	Rr 0.646 Rr 0.648 Rr 0.648	72.64 µg/ml 73.35 µg/ml 73.25 µg/ml 95.19 µg/ml 95.19 µg/ml	217.9 ng 220.0 ng 219.8 ng CV=2.20 % CV=2.20 %	(3 applications) (3 replicas)
^	Track 14 Track 15 Sample 'Glu-0413 A Volume: 3.0 µl Track 16	Rr 0.646 Rr 0.648 Rr 0.648 21-01-05' Rr 0.650	72.64 µg/ml 73.35 µg/ml 73.25 µg/ml 95.19 µg/ml 95.19 µg/ml 92.82 µg/ml	217.9 ng 220.0 ng 219.8 ng CV=2.20 % CV=2.20 % 278.5 ng	(3 applications) (3 replicas)
^	Track 14 Track 14 Sample 'Glu-0413 Volume: 3.0 µl Track 16 Track 16	Rr 0.646 Rr 0.648 Rr 0.648 Rr 0.648 Rr 0.650 Rr 0.650	72.64 µg/ml 73.35 µg/ml 73.25 µg/ml 95.19 µg/ml 95.19 µg/ml 92.82 µg/ml 96.80 µg/ml	217.9 ng 220.0 ng 219.8 ng CV=2.20 % CV=2.20 % 278.5 ng 290.4 ng	(3 applications) (3 replicas)

Spiked with 25 μ g/mL= 106.1% Spiked with 50 μ g/mL= 98.9% Spiked with 75 μ g/mL= 95.4%

Robustness







Heating of different plates at 110 °C; 120 °C and 130 °C

Calibration curves



Glc and DP4: 200 μg/mL to 5 μg/mL **DP2, DP3, DP5-DP8:** 100 μg/mL to 5 μg/mL



Analysis of the sugars in three fermenter vessels



Analysis of the sugars



Automation of the process





Parameters compared	HPLC	HPTLC
Sugars separation	Glc, DP ₂ , DP ₃ and DP ₄ ⁺	Glc, DP ₂ , DP ₃ , DP ₄ , DP ₅ , DP ₆ , DP ₇ and DP ₈ ⁺
Time for 20 samples	~600 min	~120 min
Amount of solvent per 20 samples	~360 mL	~40 mL

For the comparison an HPLC method of 30 min and a flow rate of 0.6 mL/min was considered.

Thank you!

