

Standard operating procedure for HPTLC

Version: 03	Effective Date: 26 Jul 2018
Revision history: Version 03: editorial changes Version 02: Note on use of 10x10 cm HPTLC plates added; details on derivatization procedure added (by immersion or automatic spraying) Version 01: Created 19 Dec 2012	
Related documents: ---	

Source

B. Meier & A. Spriano, Zurich University of Applied Sciences

Purpose

This standard operating procedure (SOP) provides general guidance for analysis by High-Performance Thin-Layer Chromatography HPTLC.

Definitions

HPTLC is performed on 20x10 cm HPTLC **glass plates** coated with silica gel 60 F₂₅₄. Suitable (preferably software-controlled) instruments must be employed for sample application, chromatogram development, derivatization, and documentation.

If no automatic developing chamber is available a 20x10 cm twin trough chamber may be used.*

NOTE: Record temperature and relative humidity in the laboratory

1. Preparation of plates

- 1.1 Obtain HPTLC plate silica gel 60 F₂₅₄ (20x10 cm). Record the batch number.
- 1.2 Inspect plate under UV 254 nm for any damage of the layer. If damage is detected discard plate.
- 1.3 With a soft pencil label the plate in the upper right corner with: your initials – date (dd/mm/yy) – consecutive number for the day. Example ER-23/02/10-001.
- 1.4 On the right side of the plate mark developing distance at 70 mm from lower edge of plate. NOTE: left-handed persons may label /mark the plate on left side.
- 1.5 Capture electronic images of the clean plate under UV 254 nm and under white light

2. Preparation of chamber (manual development only)

- 2.1 Obtain a twin trough chamber for 20x10 cm plates.
- 2.2 Fit the rear trough of chamber with a filter paper of corresponding size.
- 2.3 Pour 20 ml of developing solvent over the filter paper into the rear trough ensuring complete wetting.
Pour a sufficient amount of developing solvent into the front trough to have a level of 5 mm.
- 2.4 Close the lid of the chamber and allow 20 min for saturation.

* NOTE: Small HPTLC plates of 10x10 cm may be used for **manual development** only. Use a twin trough chamber for 10x10 cm plates and the following application parameters:

band length: 8 mm	sample solvent type for application: methanol
number of tracks: 7	developing solvent in rear trough (with filter paper): 10 mL
first application position X: 15 mm	developing solvent in front trough: a sufficient amount to have a level of 5 mm
application position Y: 8 mm	
distance between tracks: automatic (min. 11 mm)	

3. Sample application

3.1 Select the following application parameters on the application device:

- band length: 8 mm
- number of tracks: 15
- first application position X: 20 mm
- application position Y: 8 mm
- distance between tracks: automatic (minimum 11 mm)
- sample solvent type: methanol

3.2 Apply the volumes specified in the method for selected herbal drug.

3.3 Use track 1 to apply the system suitability test as specified in the individual method of analysis.

3.3. Disable any unused tracks.

4. Plate conditioning (manual development only)

After sample application place the plate for 45 min in a suitable desiccator containing a saturated solution of MgCl₂.

5a. Manual development

5a.1 Slowly open the lid of the saturated chamber and insert the conditioned plate into the front trough so that the back of the plate rests against the front wall of the chamber and the layer faces the inside of the chamber. Close the lid.

5a.2 Let the mobile phase ascend until it reaches the mark.

5a.3 Open the lid and remove the plate. Place it upright in a rack under a fume hood.

5a.4 Dry plate with cold air from a hair dryer for 5 min placed at a distance of 30 cm.

5b. Automatic development

Use the following settings of the automatic chamber:

- Enable pre-drying
- Saturation with filter paper 20 min
- Humidity control 10 min with MgCl₂
- Migration distance 70 mm
- Drying time 5 min
- 10 mL of developing solvent
- 25 mL of saturation solvent

NOTE: if no humidity control is available, follow step 4.

6. Derivatization

Use the reagent specified in the individual monograph.

Derivatization is either performed by using an automatic dipping device (typical immersion speed 5 cm/s and dwell time 0 sec) or an automatic derivatization device (typically 2-4 mL of reagent).

Example:

6a. Flavonoids – Derivatization by dipping

6a.1 Heat the dry plate for 5 min at 100°C.

6a.2 While hot dip plate for 1 sec into a solution of 0.5% NP reagent in ethyl acetate. Then, after two min of waiting, dip the plate for 1 sec into a solution of 5% macrogol 400 in dichloromethane.

6b. Flavonoids – Derivatization by automatic derivatization

6b.1 Heat the dry plate for 5 min at 100°C.

6b.2 After cooling to room temperature derivatize the plate with 4 mL of a 1:1 mixture of 0.5% NP reagent in ethyl acetate and 5% macrogol 400 in dichloromethane.

30 min after completion of derivatization, take an image under UV366 nm.

7. Documentation

After development, take one image each under UV 254 nm, UV 366 nm, and white (transmission + reflection). After derivatization, take one image each under UV 366 nm and white light (transmission + reflection).

8. Reporting

Create a copy of the software based report or use own reporting documents.